Abstract #2062

Pre-analytical conditions strongly influence molecular pattern in clinical biospecimen

Introduction

The development of personalized medicine in oncology is fostered by high-throughput analysis of molecular biomarkers in human cancer biospecimens. Tissue quality strongly depends on the pre-analytical conditions in which it was acquired. Therefore, warm and cold ischemia times that tissues are exposed to during collection are of critical importance. Insufficient quality of such tissues may lead to spurious results and data misinterpretation resulting in biased stratification of patients. The present study was conducted to gain a better understanding of the effects of warm and cold ischemia on the molecular composition of tissue specimens and especially on critical clinical biomarkers whose expression and activity level inform targeted therapy evaluation in cancer. NanoPro™1000 technology enables the differential quantification of protein isoforms and their phosphorylation status in a single measurement using very small amounts of protein. Protein can be detected at picogram levels, allowing analysis of post-translational modification in samples where limited protein is available (e.g., samples from small-needle biopsies) and making this a valid technology for the clinical setting. NanoPro™1000 has the potential to characterize protein phosphorylation of multiple sites using just one pan-antibody.

Methods

Patient recruitment:

Fifty patients with colorectal cancer (CRC) and 43 with hepatic metastasis of CRC who were scheduled for tumor resection surgery gave informed consent to be enrolled in the study. Only patients with a tumor larger than 3 cm in diameter were enrolled. Neoadjuvant treated patients were excluded. The study received approval by the competent ethics review committee of the medical association Hamburg.

RNA Extraction:

RNA was prepared in duplicates from every tissue block. RNA was isolated by two steps chloroform extraction phenol using followed by the RNeasy MinElute Cleanup Kit from Qiagen according to manual instruction

Gene expression:

RNA samples were analyzed in biological replicates using Affymetrix chips (whole genome; Affymetrix-Power-Tools). Statistical analysis of change in gene expression was performed by cluster analysis and direct group comparison.

Immunohistochemistry:

Tissues were formalin fixed, paraffin embedded and cut into five μm thick slices. IHC stainings were conducted using the Benchmark[®] Ultra (Roche Diagnostics Deutschland GmbH).

NanoPro™1000 Technology:

Analysis of signaling proteins was conducted using the NanoPro1000 technology platform. This technology enables the identification of multiple isoforms according to their isoelectric point. Therefore, protein lysates were separated on a nested Premix G2 5-8 gradient against the pl standard ladder 3 and immobilized for 70 seconds.

Meso Scale Discovery:

Quantification of proteins was conducted using 96-well format technology platform from Meso Scale Discovery (MSD[®], Gaithersburg, MD, USA). MSD[®] 96-well MULTI-SPOT[®] assays were performed for key proteins from the MAPK and PI3K signalling pathways.

Workflow

Study design – Sample collection

Procedure	Colonoscopy (n=48):Procedure6 biopsies/patient		Tumor resectio Removal of entire tumor and a dissection of 18 tissue pier				
Tissue type	3 x Tumor Tissue	3 x Normal Tissue		9 x Tumor Tissue			
Cold ischemia time	3 x <1 min	3 x <1 min	3 x 10 min	3 x 20 min	3 x 45 min		
Type of preservation	1 x FFPE (n=47)	1 x FFPE (n=47)	1 x FFPE (n=43)	1 x FFPE (n=43)	1 x FFPE (n=43)	_	
	2 x FF (n=48)	2 x FF (n=48)	2 x FF (n=49)	2 x FF (n=49)	2 x FF (n=49)		
		Histo	logical quality co	ntrol and sele	ction of 40 pat	ien	
Immunohisto- chemistry	1 x FFPE (n=38)	1 x FFPE (n=37)	1 x FFPE (n=39)	1 x FFPE (n=39)	1 x FFPE (n=40)	-	
Protein quantification	1 x FF (n=39)	1 x FF (n=38)	1 x FF (n=40)	1 x FF (n=40)	1 x FF (n=40)	_	
,		4	4	1			

Sample collection: etermination of tumor content normal and tumor tissue

analyses were conducted with duplicate samples



Figure 1: Flow chart showing the essential steps of tissue collection preparation and subsequent analysis. This scheme of collection and processing was applied to the colon as well as liver tissue samples. Samples for NanoPro™1000 analysis were adjusted in tumor content.

NanoPro[™]1000 - Principle of technology



Figure 2: The NanoPro 1000 Technology enables the identification of multiple isoform phosphorylations by separating low amounts of proteins according to their isoelectric point. Therefore, by integrating this technology into our drug testing platform we are able obtain more detailed information about activation and inactivation of signaling proteins. ([™] ProteinSimple)

Unger, F.T.; Lange, N.; Uhlig, P.C.; Krüger, J.; Juhl, H. and David, K.A.

Indivumed GmbH, Hamburg, Germany

This project was funded in part with federal funds from the National Cancer Institute, National Institutes of Health, under contract HHSN261200800001E

Normal colon tissue



Changes in gene expression due to ischemia (Affymetrix)

			Pre vs 10'		Pre v	
			Log-fold			
	Gene	Protein and function	p-value	change	p-value	
	CYR61	Cysteine-rich angiogenic inducer 61; extracellular matrix- associated signaling protein that plays important roles in tissue repair	1.57E-18	2.49	2.84E-22	
	RGS1	Regulator of G-protein signaling 1; attenuates signaling activity	4.21E-21	2.34	5.05E-25	
		of G-proteins	1.65E-22	2.33	8.96E-25	
	EGRI	Early growth response 1; transcription factor	1.92E-12	2.02	1.79E-18	
		-	1.4E-10	1.98	2.57E-16	
	SGK1	Serum/glucocorticoid regulated kinase 1; activates potassium, sodium and chloride channels	4.03E-14	2.09	2.74E-16	
	FOS	FBJ murine osteosarcoma viral oncogene homolog; transcription factor involved in cell proliferation, differentiation, survival, hypoxia and angiogenesis	1.26E-18	2.77	9.3E-25	
	C8orf4	Chromosome 8 open reading frame 4; uncharacterized protein	2.39E-12	2	3.64E-15	
	DUSP1	Dual specificity phosphatase 1; dephosphorylates MAP kinase MAPK1/ERK2	1.21E-17	1.85	3E-23	
	DUOX2	Dual oxidase 2; plays a role in antimicrobial defense at the mucosal surface	3.36E-08	- <mark>1.88</mark>	2.24E-10	
	SLC6A14	Solute carrier family 6 (amino acid transporter), member 14; mediates the uptake of a broad range of amino acids	0.0000101	-1.87	0.0000036	
	VNN1	Vanin 1; amidohydrolase recycling pantothenic acid (vitamin B5) and releasing cysteamine	0.000102	-1.62	0.000000164	

Table 1: Differentially expressed genes in normal and colorectal tumor tissue. Gene expression was compared: pre, before hepatic pedicle clamping; post, after clamping; 10[°], 10 minutes after resection; 20[°], 20 minutes after resection, and 45[°], 45 minutes after resection.



Figure 5: Time-dependent changes in overall phosphorylation of AKT (A), ERK1/2 (B) MEK1/2 (C), and c-MET (D) in normal (N) compared with tumor (T) Tissue from 20 patients are shown by box-and-whisker plots. Results were displayed as fold changes normalized in normal tissue to sample NC4/NO and in tumor tissue to sample NC4/TO Kruskal–Wallis test and Dunn test for multiple comparisons or analysis of variance (ANOVA) were used for statistical analysis. N = normal tissue; T = tumor tissue; 0 =before surgery; 10, 20, 45 = 10, 20, 45 min after resection; $a=p\leq0.05$; $b=p\leq0.01$.







Figure 6: Changes in isoform phosphorylation of the target signaling proteins, AKT, ERK1/2, MEK1/2 and c-MET in response to ischemia. An overview of all isoforms of AKT. . ERK1/2. MEK1/2 and c-MET in normal tissue (A-D) of all patients analyzed is shown: there results were displayed as fold changes normalized in normal tissue to sample NC5/NO and in tumor tissue (E-H) to sample NC5/TO. An overview of all isoforms in tumor tissue of all patients analyzed is shown in below. Kruskal–Wallis test and Dunn test for multiple comparisons or analysis of variance (ANOVA) were used for statistical analysis. N = normal tissue; T = tumor tissue; T

